

# Triple-bagged Pipets for Clean Room Applications

## Falcon® and Stripette® Serological Pipets




Triple-bagged packaging

### Triple-bagged for Safety and Convenience

Corning's triple-bagged pipets offer the safest and most convenient solution for applications in which sterility needs to be assured. This packaging not only protects the product from external contamination, but also features a Sterility Assurance Level (SAL) of  $10^{-6}$  and is suitable for laboratory clean rooms that require multiple packaging layers.

- ▶ Triple-bagged
- ▶ Each pipet is individually wrapped in easy open plastic/plastic or paper/plastic
- ▶ No need for decontamination spraying with alcohol
- ▶ RNase-/DNase-free, and human DNA-free
- ▶ Sterility Assurance Level (SAL) of  $10^{-6}$
- ▶ Extractables data available upon request.

**NOTE:** Triple-bagged Falcon and Stripette pipets have the same features as standard pipets except for the packaging.



**1**  
Contained in a twist-tied liner (Bag 1), the entire case contents of triple-bagged pipets can easily be removed from the corrugated box for clean area storage.

**2**  
At the first stage, the liner (Bag 1) is opened and the double-bagged inner packs with their outer sterile barriers (Bag 2) are removed.

**3**  
At the second stage the outer sterile barrier (Bag 2) is opened and the inner sterile barrier (Bag 3) is removed.

**4**  
Reaching the point of use, the inner sterile barrier (Bag 3) is opened and the individually wrapped pipets are removed.

### Ordering Information

#### Falcon Triple-bag Packaging, Individually Wrapped, Plastic/Plastic\*

Cat. No.	Capacity (mL)	Graduations (mL)	Negative Grads. (mL)	Color-coded Stripe	Qty/Pk	Qty/Cs
355024*	1	1/100	0.2	Yellow	10	1,000
355025*	2	1/100	0.2	Green	10	1,000
355026*	5	1/10	2.5	Blue	10	200
355027*	10	1/10	3.0	Red	10	200
355028*	25	2/10	10.0	Purple	10	200
355029*	50	1/2	10.0	Black	10	100

\*Minimum Order Quantity. Please contact your local Corning Account Manager or Corning Customer Service for lead time and minimum order quantity required.

#### Stripette Triple-bag Packaging, Individually Wrapped, Paper/Plastic

Cat. No.	Capacity (mL)	Graduations (mL)	Negative Grads. (mL)	Color-coded Stripe	Qty/Pk	Qty/Cs
7041	1	1/100	0.2	Yellow	50	1,000
7042	2	1/100	0.2	Green	50	1,000
7045	5	1/10	2.5	Blue	50	200
7015	10	1/10	3.0	Orange	50	200
7016	25	2/10	10.0	Red	25	200
7017	50	1/2	10.0	Purple	25	100
7000	100	1	N/A	Aqua	10	100

**Warranty/Disclaimer:** Unless otherwise specified, all products are for research use or general laboratory use only.\* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. \*For a listing of US medical devices, regulatory classifications or specific information on claims, visit [www.corning.com/resources](http://www.corning.com/resources).

# CORNING

Corning Incorporated  
Life Sciences

[www.corning.com/lifesciences](http://www.corning.com/lifesciences)

## **NORTH AMERICA**

t 800.492.1110  
t 978.442.2200

## **ASIA/PACIFIC**

### **Australia/New Zealand**

t 61 427286832

### **Chinese Mainland**

t 86 21 3338 4338

### **India**

t 91 124 4604000

### **Japan**

t 81 3-3586 1996

### **Korea**

t 82 2-796-9500

### **Singapore**

t 65 6572-9740

### **Taiwan**

t 886 2-2716-0338

## **EUROPE**

CSEurope@corning.com

### **France**

t 0800 916 882

### **Germany**

t 0800 101 1153

### **The Netherlands**

t 020 655 79 28

### **United Kingdom**

t 0800 376 8660

## **All Other European Countries**

t +31 (0) 206 59 60 51

## **LATIN AMERICA**

grupolA@corning.com

### **Brazil**

t 55 (11) 3089-7400

### **Mexico**

t (52-81) 8158-8400